



VirClia® Lotus

The revolution in infectious serology



vircell 
MICROBIOLOGISTS

WIDEST PANEL FOR INFECTIOUS SEROLOGY IN CLIA



BACTERIAL RESPIRATORY INFECTIONS

- VCM092 BORDETELLA PERTUSSIS TOXIN VIRCLIA® IgG MONOTEST (quant.)
VCM016 CHLAMYDOPHILA PNEUMONIAE VIRCLIA® IgA MONOTEST
VCM017 CHLAMYDOPHILA PNEUMONIAE VIRCLIA® IgG MONOTEST
VCM018 CHLAMYDOPHILA PNEUMONIAE VIRCLIA® IgM MONOTEST
VCM019 COXIELLA BURNETII VIRCLIA® IgG MONOTEST
VCM020 COXIELLA BURNETII VIRCLIA® IgM MONOTEST
VCM049 LEGIONELLA PNEUMOPHILA SG 1 VIRCLIA® IgG MONOTEST
VCM050 LEGIONELLA PNEUMOPHILA SG 1 VIRCLIA® IgM MONOTEST
VCM052 LEGIONELLA PNEUMOPHILA SG 1-6 VIRCLIA® IgG+IgM MONOTEST
VCM062 MYCOPLASMA PNEUMONIAE VIRCLIA® IgG MONOTEST
VCM064 MYCOPLASMA PNEUMONIAE VIRCLIA® IgM MONOTEST



VIRAL RESPIRATORY INFECTIONS

- VCM002 ADENOVIRUS VIRCLIA® IgG MONOTEST
VCM097 COVID-19 VIRCLIA® IgG MONOTEST **NEW**
VCM098 COVID-19 VIRCLIA® IgM+IgA MONOTEST **NEW**
VCM100 COVID-19 SPIKE QUANTITATIVE VIRCLIA® IgG MONOTEST **NEW**
VCM042 INFLUENZA A VIRCLIA® IgG MONOTEST
VCM045 INFLUENZA B VIRCLIA® IgG MONOTEST
VCM065 PARAINFLUENZA 1 VIRCLIA® IgG MONOTEST
VCM068 PARAINFLUENZA 2 VIRCLIA® IgG MONOTEST
VCM071 PARAINFLUENZA 3 VIRCLIA® IgG MONOTEST
VCM077 RESPIRATORY SYNCYTIAL VIRUS VIRCLIA® IgG MONOTEST



VACCINE-PREVENTABLE INFECTIONS

- VCM006 DIPHTHERIA VIRCLIA® IgG MONOTEST (quant.)
VCM054 MEASLES VIRCLIA® IgG MONOTEST
VCM057 MEASLES VIRCLIA® IgM MONOTEST
VCM060 MUMPS VIRCLIA® IgG MONOTEST
VCM061 MUMPS VIRCLIA® IgM MONOTEST
VCM004 TETANUS VIRCLIA® IgG MONOTEST (quant.)
VCM089 VARICELLA-ZOSTER VIRCLIA® IgG MONOTEST
VCM091 VARICELLA-ZOSTER VIRCLIA® IgM MONOTEST



VECTOR-BORNE INFECTIONS

- VCM009 BORRELIA VIRCLIA® IgG MONOTEST
VCM010 BORRELIA VIRCLIA® IgM MONOTEST
VCM008 CHAGAS VIRCLIA® IgG+IgM MONOTEST
VCM099 CHAGAS TESA VIRCLIA® IgG+IgM MONOTEST **NEW**
VCM059 CHIKUNGUNYA VIRCLIA® IgG MONOTEST
VCM063 CHIKUNGUNYA VIRCLIA® IgM MONOTEST
VCM023 DENGUE VIRCLIA® IgG MONOTEST
VCM024 DENGUE VIRCLIA® IgM MONOTEST (capture)
VCM053 LEISHMANIA VIRCLIA® IgG+IgM MONOTEST
VCM080 RICKETTSIA CONORII VIRCLIA® IgG MONOTEST
VCM082 RICKETTSIA CONORII VIRCLIA® IgM MONOTEST
VCM069 ROSS RIVER VIRUS VIRCLIA® IgG MONOTEST
VCM070 ROSS RIVER VIRUS VIRCLIA® IgM MONOTEST
VCM047 TICK-BORNE ENCEPHALITIS VIRCLIA® IgG MONOTEST
VCM048 TICK-BORNE ENCEPHALITIS VIRCLIA® IgM MONOTEST
VCM078 WEST NILE VIRUS VIRCLIA® IgG MONOTEST
VCM079 WEST NILE VIRUS VIRCLIA® IgM MONOTEST
VCM033 ZIKA VIRCLIA® IgG MONOTEST
VCM035 ZIKA VIRCLIA® IgM MONOTEST



FUNGAL INFECTIONS

- VCM073 ASPERGILLUS GALACTOMANNAN Ag VIRCLIA® MONOTEST
VCM094 INVASIVE CANDIDIASIS (CAGTA) VIRCLIA® IgG MONOTEST



SEXUALLY TRANSMITTED INFECTIONS

- VCM012 CHLAMYDIA TRACHOMATIS VIRCLIA® IgA MONOTEST
VCM013 CHLAMYDIA TRACHOMATIS VIRCLIA® IgG MONOTEST
VCM015 CHLAMYDIA TRACHOMATIS VIRCLIA® IgM MONOTEST
VCM085 SYPHILIS VIRCLIA® IgG MONOTEST
VCM086 SYPHILIS VIRCLIA® IgG+IgM MONOTEST
VCM032 SYPHILIS VIRCLIA® IgM MONOTEST



TORCH INFECTIONS

- VCM038 CYTOMEGALOVIRUS VIRCLIA® IgG AVIDITY MONOTEST
VCM021 CYTOMEGALOVIRUS VIRCLIA® IgG MONOTEST (quant.)
VCM022 CYTOMEGALOVIRUS VIRCLIA® IgM MONOTEST (capture)
VCM031 HERPES SIMPLEX 1 VIRCLIA® IgG MONOTEST
VCM039 HERPES SIMPLEX 1 VIRCLIA® IgM MONOTEST
VCM034 HERPES SIMPLEX 1+2 VIRCLIA® IgG MONOTEST
VCM036 HERPES SIMPLEX 1+2 VIRCLIA® IgM MONOTEST
VCM037 HERPES SIMPLEX 2 VIRCLIA® IgG MONOTEST
VCM041 HERPES SIMPLEX 2 VIRCLIA® IgM MONOTEST
VCM083 RUBELLA VIRCLIA® IgG MONOTEST (quant.)
VCM084 RUBELLA VIRCLIA® IgM MONOTEST (capture)
VCM025 TOXOPLASMA VIRCLIA® IgG AVIDITY MONOTEST
VCM087 TOXOPLASMA VIRCLIA® IgG MONOTEST (quant.)
VCM088 TOXOPLASMA VIRCLIA® IgM MONOTEST (capture)



ZOONOTIC INFECTIONS

- VCM046 BARTONELLA HENSELAE VIRCLIA® IgG MONOTEST
VCM043 BARTONELLA HENSELAE VIRCLIA® IgM MONOTEST
VCM005 BRUCELLA VIRCLIA® IgG MONOTEST
VCM007 BRUCELLA VIRCLIA® IgM MONOTEST
VCM040 HYDATIDOSIS VIRCLIA® IgG MONOTEST
VCM011 TULAREMIA VIRCLIA® IgG+IgM MONOTEST



OTHER INFECTIONS

- VCM044 EPSTEIN-BARR EA VIRCLIA® IgG MONOTEST
VCM026 EPSTEIN-BARR EBNA VIRCLIA® IgG MONOTEST
VCM027 EPSTEIN-BARR VCA VIRCLIA® IgG MONOTEST
VCM028 EPSTEIN-BARR VCA VIRCLIA® IgM MONOTEST
VCM029 HELICOBACTER PYLORI VIRCLIA® IgA MONOTEST
VCM030 HELICOBACTER PYLORI VIRCLIA® IgG MONOTEST
VCM066 HEPATITIS E VIRCLIA® IgG MONOTEST (quant.)
VCM067 HEPATITIS E VIRCLIA® IgM MONOTEST (capture)
VCM003 HHV-6 VIRCLIA® IgG MONOTEST
VCM001 INFECTIOUS MONONUCLEOSIS HETEROPHILE VIRCLIA® IgM MONOTEST
VCM014 LEPTOSPIRA VIRCLIA® IgM MONOTEST
VCM075 PARVOVIRUS VIRCLIA® IgG MONOTEST
VCM076 PARVOVIRUS VIRCLIA® IgM MONOTEST (capture)



VirClia® Lotus

Any test, any sample, any time

RANDOM ACCESS, CONTINUOUS LOADING

VirClia® Lotus accepts any combination of samples and reagents, regardless of the instrument status. Even while running, it remains accessible for reagent and sample addition.

STAT FUNCTION

Urgent samples can be prioritized. With VirClia® Lotus, STAT test requests are given the highest priority for processing.

TOTAL FLEXIBILITY

This true monotest solution allows the laboratory to run esoteric and low-volume infectious serology tests daily, alongside the laboratory routine. VirClia® Lotus unifies in just one CLIA system, all those parameters typically performed with ELISA, IFA or Rapid Tests.

Positive identification of samples



QR code ID of the strips



High precision single tip



POSITIVE IDENTIFICATION OF SAMPLES & STRIPS

VirClia® Lotus built-in scanner allows for full traceability of both samples and strips, avoiding clerical errors.

SAMPLE PRIMARY TUBES

Save hands-on-time using primary tubes (11-16 mm) directly in VirClia® Lotus, avoiding the inconvenience of transferring samples or performing external predilutions.

EXCELLENT TURNAROUND TIME

Results appear as soon as they are processed. The first results are available within 70 minutes and then, every 35 seconds thereafter.

The system is able to process 40 tests in 90 minutes. VirClia® Lotus holds up to 79 different monotests on-board and 50 primary sample tubes.

VIRCLIA® THE TRUE MONOTEST

All In One

Vircell provides a complete, fully automated platform together with self-contained monotests assays with all the needed reagents in the strip including calibrators and controls.

All Are One

About 90 tests available with compatible protocols, no special workflow considerations. All different assays can be treated as one.



NO HIDDEN COSTS

- > On-demand testing, 1 monotest = 1 reportable test
- > Individual quality control per monotest, no need for extra controls or calibrations
- > Shelf life associated to each individual strip; no throwing away unused components after opening!

TRUE SELF-CONTAINED MONOTEST

- All reagents are included and ready-to-use
- Every monotest includes a calibrator for internal validation
- QR code ID of the strip
- One strip= one result. Nothing else is required

MULTIPARAMETER

- Unmatched infectious diseases menu in CLIA
- Allows to run any number of different VirClia® parameters simultaneously
- Continuous loading, multi-batch testing
- STAT testing capability



CLIA, CHEMILUMINESCENT IMMUNOASSAY

- Highly sensitive detection method
- Extended linearity range
- Qualitative and quantitative tests
- Excellent turnaround time

VIRCLIA® LOTUS TECHNICAL DATA

> OVERVIEW

Method	Chemiluminescence immunoassay (CLIA)
Operation mode	Continuous loading, random access
Urgent Samples Management	STAT function
Multiparameter	Full compatibility with all VirClia® assays

> CAPACITY

Sample loading	50 Primary tubes, 50 dilution tubes and 6 bottles of diluent
Strip loading	39 in loading rotor and 40 in reaction rotor
Wash buffer	5 liters
Distilled water	0.25 liters
Liquid waste	5 liters (internal tank) / 11 liters (external tank)
Solid waste	80 strips (automatic disposal)

> THROUGHPUT

Time to first result	First results in 70 min (one result every 35 seconds)
Speed	40 tests in 90 min

> SAMPLES AND REAGENTS HANDLING

Sample tubes	Primary tube (11-16 mm diameter)
Predilution	Automatic predilution if required
Built-in scanner	Positive identification of samples and strips

> USER INTERFACE

PC	External laptop 15.3" Touch screen, 1920 x 1080 resolution
Operating system	Windows 10 Pro x64
LIS connectivity	Bidirectional VirCom® middleware
Results management	Customized mode (per test, per patient) Automatic full report (QC per test, no standard curves) Accessible results

> OTHER FEATURES

Dimensions and weight	W 75 cm x D 74 cm x H 65 cm. 70 kg (max)
Noise Level	<75 decibels
Power supply	Europe: 230 Vac@ 50Hz USA/Canada: 110-120Vac@ 60Hz
Electrical Power Consumption	365 VA

